

sIPV shipping box as per *WHO Guidelines on the International Packaging and Shipping of Vaccines, Sixth Edition*

<b>Type</b>	<b>Dimensions (mm) WxLxLH</b>	<b>Loading method</b>	<b>Number of vials</b>	<b>Total weight (kg)</b>
LCX-96VP	660 x 490 x 600	/	1,800	49.95
XT12-930PH	1100 x 995 x 1560	Single-layer coolant	7,200	289
		Double-layer coolant	6,400	348.4
YT-TP-6 (1 vial/ carton)	1100 x 1000 x 1500	Single-layer coolant	7,200	273.8
		Double-layer coolant	5,600	344.6
YT-TP-6 (3 vials/ carton)	1100 x 1000 x 1500	Single-layer coolant	20,400	355.34
		Double-layer coolant	16,200	410.97
YT-VIP-4 (1 vial/ carton)	850 x 730 x 670	Single-layer coolant	1,200	84.1
		Double-layer coolant	800	104.6
YT-VIP-4 (3 vials/ carton)	850 x 730 x 670	Single-layer coolant	3,600	99.46
		Double-layer coolant	2,400	114.84