(Desano Pharmaceuticals Private Ltd), HA746

This part outlines the scientific assessment and knowledge about this product at the time of prequalification. Updates to this information are included in parts 1 to 5 and 8 of this WHOPAR.

## **SCIENTIFIC DISCUSSION**

Name of the Finished Pharmaceutical Product	[HA746 trade name]*
Manufacturer of Prequalified Product	Shanghai Desano Bio-Pharmaceutical Co., Ltd., 1479, Zhangheng Road, China (Shanghai) Pilot Free Trade Zone, Shanghai, 201203, P.R. China
Active Pharmaceutical Ingredient(s) (API)	Dolutegravir (as sodium)/lamivudine/tenofovir disoproxil fumarate
Pharmaco-therapeutic group (ATC Code)	Antivirals for treatment of HIV infections, combinations, (J05AR27)
Therapeutic indication	[HA746 trade name] is indicated for the treatment of human immunodeficiency virus (HIV) infection in adults and adolescents weighing at least 30 kg.

### 1. Introduction

[HA746 trade name] is indicated for the treatment of human immunodeficiency virus (HIV) infection in adults and adolescents weighing at least 30 kg. Detailed information on the use of this product is described in the summary of product characteristics (SmPC).

[HA746 trade name] should be initiated by a health care provider experienced in the management of HIV infection.

## 2. Assessment of quality

The assessment was done in accordance with the requirements of WHO's Guidelines on submission of documentation for a multisource (generic) finished pharmaceutical product for the WHO Prequalification of Medicines Programme: quality part.

# **Active pharmaceutical Ingredient (API)**

Dolutegravir sodium, lamivudine and tenofovir disoproxil fumarate (TDF) have been prequalified by WHO according to WHO's *Procedure for assessing the acceptability, in principle, of active pharmaceutical ingredients for use in pharmaceutical products* (WHO Technical Report Series No. 953, 2009, Annex 4). This procedure provides an assurance that the APIs, used in the manufacture of [HA746 trade name], are of good quality and manufactured in accordance with WHO Good Manufacturing Practices. API prequalification consists of a comprehensive evaluation procedure that has two components: Assessment of the API master file (APIMF) to verify compliance with WHO norms and standards, and assessment of the sites of API manufacture to verify compliance with WHO GMP requirements.

<sup>\*</sup> Trade names are not prequalified by WHO. This is the national medicines regulatory authority's responsibility.

Page 1 of 5

Dolutegravir (sodium)/lamivudine/tenofovir disoproxil fumarate 50mg/300mg/300mg tablets

(Desano Pharmaceuticals Private Ltd), HA746

#### Other ingredients

Other ingredients used in the core tablet formulation include lactose monohydrate, microcrystalline cellulose, pregelatinized starch, croscarmellose sodium, magnesium stearate, mannitol, sodium starch glycolate, povidone and pigment yellow 42. The commercially sourced proprietary pigment yellow which is included in the tablet formulation is supported by appropriate declarations and controlled by acceptable specifications. The commercially sourced proprietary film-coating mixture contains polyvinyl alcohol-partially hydrolysed, talc, macrogol/PEG and titanium dioxide. TSE / BSE free certificates have been provided for the excipients.

## Finished pharmaceutical product (FPP)

## Pharmaceutical development and manufacture

The multisource product is a white coloured, capsule-shaped tablet with 'D15' debossed on lower side and plain on the upper side. The tablet should be free from physical defects. The tablets are packaged in a white, round, opaque HDPE bottle with white, polypropylene child resistant closure with white outer shell embossed with pictorial design and clear inner shell fitted with a heat seal and pulp liner. The bottle also contains either one or two silica gel sachets as desiccant.

The objective of the development programme was to obtain a stable, robust, immediate-release FDC tablet that is bioequivalent to the WHO recommended comparator products: Tivicay® (dolutegravir) 50mg tablets, Epivir® (lamivudine) 300mg tablets and Viread® (TDF) 300mg tablets. The selection of excipients was based on their suitability to achieve the desired tablet characteristics, information of the qualitative composition of the comparator products, compatibility with the APIs and literature studies. A bilayer tablet containing the BCS low soluble dolutegravir sodium in one layer and the highly soluble lamivudine and TDF in the other layer was developed. Wet granulation technique was selected as the most ideal process based on the high content and poor compressibility of the APIs. Formulation trials were performed to optimise the concentration of excipients and process parameters. Satisfactory in-process controls have been established.

According to a risk evaluation by the applicant, the FPP appears to have no potential to contain nitrosamine impurities and hence no risk was identified.

## **Specifications**

The finished product specifications include tests for description, identification of the APIs (HPLC and UV with PDA detector) and colorants, average weight, water content (KF), uniformity of dosage units (by content uniformity), dissolution (HPLC detection), assay (HPLC), related substances (HPLC), residual solvents (GC), elemental impurities and microbial limits. The test methods have been satisfactorily validated.

#### Stability testing

Stability studies have been conducted at 30°C/75%RH as long-term storage conditions and for six months at accelerated conditions in the packaging proposed for marketing of the product. The data showed slight degradation for TDF, though all parameters were well within the agreed limits at both storage conditions. Based on the available stability data, the proposed shelf life and storage conditions as stated in the SmPC are acceptable. The in-use storage period after first opening of the bottle (90's HDPE container) is based on in-use stability data.

#### Conclusion

The quality part of the dossier is accepted.

(Desano Pharmaceuticals Private Ltd), HA746

#### 3. Assessment of bioequivalence

The following bioequivalence study has been performed in 2019 according to internationally accepted guidelines:

A randomized, open-label, balanced, two-treatment, two-period, two-sequence, single-dose, crossover oral bioequivalence study of test product [HA746 trade name] of Aizant Drug Research Solutions Pvt. Ltd., India with Reference product (R= R1: EPIVIR® (lamivudine) 300 mg tablets manufactured by GlaxoSmithKline for ViiV Healthcare (ViiV)+ R2: VIREAD® (tenofovir disoproxil fumarate) 300 mg tablets manufactured by Gilead Sciences, Inc. (Gilead) + R3: TIVICAY® (dolutegravir) 50 mg tablets manufactured by GlaxoSmithKline for ViiV Healthcare (ViiV) in normal healthy adult human subjects under fasting conditions (study no. C17276).

The objective of the study was to compare the bioavailability of the stated Dolutegravir/Lamivudine/Tenofovir Disoproxil Fumarate 50mg/300mg/300mg tablet manufactured for Aizant Drug Research Solutions Pvt. Ltd., India (test drug) with the reference formulations Tivicay® (ViiV Healthcare), Epivir® (ViiV Healthcare) and Viread® (Gilead Sciences, Inc.) and to assess bioequivalence. The comparison was performed as a single centre, open label, randomized, crossover study in healthy subjects under fasting conditions. Each subject was assigned to receive each of the following two treatments in a randomized fashion:

Treatment T: Test -1 tablet [HA746 trade name]

(dolutegravir 50 mg + lamivudine 300 mg + tenofovir disoproxil fumarate 300

mg)

Batch no. DLTTA1802B.

Treatment R: Reference

− 1 tablet Tivicay® (dolutegravir 50 mg)

Batch no. 7ZP6705

− 1 tablet Epivir® (lamivudine 300 mg)

Batch no. 8ZP8167.

− 1 tablet Viread<sup>®</sup> (tenofovir disoproxil fumarate 300 mg)

Batch no. 007674.

A minimal 12-day wash-out period was observed between administration of test and references. Serial blood samples (1 pre-dose sample and 24 samples within 72h post dose) were taken during each study period to obtain bioavailability characteristics AUC,  $C_{max}$  and  $t_{max}$  for bioequivalence evaluation. Drug concentrations for dolutegravir, lamivudine and tenofovir were analyzed using validated LC-MS/MS methods. The limit of quantification was stated to be about 20 ng/mL for dolutegravir, 15 ng/mL for lamivudine and 3 ng/mL for tenofovir.

The study was performed with 62 participants; data generated from a total of 57 subjects were utilized for analysis to establish pharmacokinetic parameters and assess bioequivalence.

Arithmetic mean and geometric mean values of the pharmacokinetic variables for dolutegravir, lamivudine and tenofovir as well as statistical results are summarised in the following tables:

#### **Dolutegravir**

	Test formulation (T)	Reference (R)	log-transformed parameters	
Pharmacokinetic	arithmetic mean ± SD	arithmetic mean $\pm$ SD	Ratio	Conventional
Parameter	(geometric mean)	(geometric mean)	T/R (%)	90% CI
				(ANOVAlog)
t <sub>max</sub> (h)	$3.03 \pm 1.84$	$2.58 \pm 1.77$	-	-
C <sub>max</sub> (ng/mL)	$3097 \pm 655$	$3027 \pm 965$	104.9	99.4 – 110.6
	(3030)	(2889)		
AUC <sub>0-t</sub> (ng·h/mL)	$65708 \pm 18381$	$60984 \pm 19836$	109.0	103.1 - 115.2
	(63376)	(58148)		
AUC <sub>0-inf</sub> (ng·h/mL)	$70057 \pm 20743$	$66047 \pm 24522$	108.2	102.0 - 114.7
	(67289)	(62214)		

(Desano Pharmaceuticals Private Ltd), HA746

#### Lamivudine

	Test formulation (T)	Reference (R)	log-transformed parameters	
Pharmacokinetic	arithmetic mean ± SD	arithmetic mean ± SD	Ratio	Conventional
Parameter	(geometric mean)	(geometric mean)	T/R (%)	90% CI
				(ANOVAlog)
t <sub>max</sub> (h)	$2.24 \pm 1.10$	$1.77 \pm 0.98$	-	-
C <sub>max</sub> (ng/mL)	$2062 \pm 452$	$2198 \pm 561$	94.8	90.0 – 99.8
	(2014)	(2126)		
AUC <sub>0-t</sub> (ng·h/mL)	$13226 \pm 2839$	$13113 \pm 2789$	100.9	96.4 – 105.7
	(12937)	(12817)		
AUC <sub>0-inf</sub> (ng·h/mL)	$13435 \pm 2826$	$13330 \pm 2776$	100.9	96.5 – 105.4
	(13153)	(13042)		

#### **Tenofovir**

	Test formulation (T)	Reference (R)	log-transformed parameters	
Pharmacokinetic	arithmetic mean $\pm$ SD	arithmetic mean $\pm$ SD	Ratio	Conventional
Parameter	(geometric mean)	(geometric mean)	T/R (%)	90% CI
				(ANOVAlog)
$t_{max}(h)$	$1.29 \pm 0.76$	$0.90 \pm 0.43$	-	-
C <sub>max</sub> (ng/mL)	$402 \pm 138$	$437 \pm 150$	92.3	85.9 – 99.1
	(381)	(413)		
$AUC_{0-t} (ng \cdot h/mL)$	$3099 \pm 734$	$3012 \pm 748$	103.4	99.4 – 107.5
	(3025)	(2926)		
AUC <sub>0-inf</sub> (ng·h/mL)	$3318 \pm 803$	$3252 \pm 799$	102.6	98.8 – 106.5
	(3238)	(3157)		

The results of the study show that preset acceptance limits of 80 -125 % are met by both AUC and  $C_{max}$ values regarding dolutegravir, lamivudine and tenofovir. Accordingly, the test [HA746 trade name] meets the criteria for bioequivalence with regard to the rate and extent of absorption and is therefore bioequivalent to the reference formulations Tivicay® (ViiV Healthcare), Epivir® (ViiV Healthcare) and Viread® (Gilead Sciences, Inc.).

## Summary of product safety and efficacy

[HA746 trade name] has been shown to conform to the same relevant standards of quality, efficacy and safety as those required of the comparator product. According to the submitted data on quality and bioavailability, [HA746 trade name] is pharmaceutically and therapeutically equivalent and thus interchangeable with the comparator product Tivicay® (ViiV Healthcare), Epivir® (ViiV Healthcare) and Viread® (Gilead Sciences, Inc.) for which benefits have been proven in terms of clinical efficacy. The clinical safety of [HA746 trade name] is considered acceptable when guidance and restrictions stated in the summary of product characteristics (SmPC) are considered. Refer to the SmPC (WHOPAR part 4) for data on clinical safety.

## 5. Benefit risk assessment and overall conclusion

## Quality

Physicochemical and biological aspects relevant to the uniform pharmaceutical characteristics have been investigated and are controlled in a satisfactory way. The quality of this product is considered to lead to an acceptable clinical performance when [HA746 trade name] is used in accordance with the SmPC.

Dolutegravir (sodium)/lamivudine/tenofovir disoproxil fumarate 50mg/300mg/300mg tablets

(Desano Pharmaceuticals Private Ltd), HA746

## Bioequivalence

[HA746 trade name] has been shown to be bioequivalent with Tivicay® (ViiV Healthcare), Epivir® (ViiV Healthcare) and Viread® (Gilead Sciences, Inc.).

## **Efficacy and Safety**

Regarding clinical efficacy and safety, [HA746 trade name] is considered effective and safe to use when the guidance and restrictions in the SmPC are taken into consideration.

#### **Benefit Risk Assessment**

Based on WHO's assessment of data on quality, bioequivalence, safety and efficacy the team of assessors considered that the benefit—risk profile of [HA746 trade name] was acceptable for the following indication: 'treatment of human immunodeficiency virus (HIV) infection in adults and adolescents weighing at least 30 kg ', and would allow inclusion of [HA746 trade name], manufactured at Shanghai Desano Bio-Pharmaceutical Co., Ltd., 1479, Zhangheng Road, China (Shanghai) Pilot Free Trade Zone, Shanghai, 201203, P.R. China in the list of prequalified medicinal products.